QiS

MANUAL

HANDLEIDING

MODE D'EMPLOI

ANLEITUNG

P822
P862

April 2006

This instrument is manufactured with the latest technology and needs no particular maintenance. Qis certifies that this instrument was thoroughly inspected and tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may differ.

The normal operating temperature should be between 0° and 40°C. Never use the instrument in a room with high humidity (>95 %) or at very low temperatures (condensation water!).

Warranty

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of thirty-six (36) months from the date of shipment ex factory. Qis will repair all defective equipment returned to it during the warranty period without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. A return authorisation must be obtained from Qis before returning any product for warranty repair on a freight prepaid basis!

Qis is not liable for consequential damages arising out of the use or handling of its products.

Servicing

In the event of this instrument being returned for servicing, the owner is requested to remove the power supply lead and **NOT** to send the following items unless they are suspect:

Manual

Cables

Accessories

If serious malfunctioning occurs, stop using the unit immediately and consult your local Qis dealer.

Keyboard


MODE = Selects all modes or escapes from error traps, calibration procedures, etc..by returning to the original mode.

←→↑↓ = Button for entering a value or for selecting a function.

CAL = Starts or proceeds a calibration or a function.

STOP = Stops the control without interrupting the measurements.

HELP = Built-in manual.

pH/mV/O₂ electrode

Both channels are provided with floating inputs. The maximum common mode voltage (potential between the (-) terminals of CHAN.1 and CHAN.2) should not exceed ±1 V. Connect the internal wire of the electrode cable to terminal (+) and the screening of the cable to terminal (-) of (CHAN1 or CHAN2). Ground loops can be avoided by immersing an earthing rod in the process liquid close to the electrodes. Connect it only to terminal (M)!

Temperature probe

A separate automatic temperature compensation is provided for each channel but a common Pt1000 (or Pt100) thermocompensator for both channels can also be used. Connect the internal wire of the electrode cable to terminal (+) and the screening of the cable to terminal (-) of (C1 or C2).

Recorder

Two recorders, one for each channel, can be connected to the (+) and (-) terminals of (OUT1) and (OUT2). Max. impedance = 300 Ω.

Relays

There are three (P862) or four (P822) relays (REL1...REL4) with voltage-free contacts. Connect their two contacts in series with the load. The ratings of the connected load should lie between 12...250 V~, maximum 0.5 A.

Mains


Connect phase to terminal (P), neutral to terminal (N), and earth to earthing terminal.

- Always short circuit the unused pH/mV/O₂ input terminals!
- Special care should be taken when connecting a typical pH/mV/O₂ coaxial cable. Make sure that all parts are dry and clean and do not forget to remove the black antistatic foil between the inner lead and the outer screening!

RS485

Only use twisted pair UTP cable to connect up to 31 controllers to a computer!

Before the exact location of the electrodes is finalised, it is worthwhile to carry out a measuring profile of the process stream to ensure that the required homogeneity is being achieved. As a general rule the electrodes should protrude into the process stream by one-third of the pipe width and be positioned so as to avoid gas pockets or sediments.

Generally the measuring electrodes are inserted into the flow with the sample pick-up facing downstream in order to reduce the risk of blockage by suspended material.

All sensors are usually positioned vertically at a suitable point downstream of the dosing point, where the measurements are stable and reliable. The response of the dosing system will be affected by the time required to disperse the reagents and their arrival at the measuring point.

Language

1. Press **SET**.
2. Select [**LANGUAGE**] and press **CAL**.
3. Select the desired language and press **CAL**.

Date and Time

1. Press **SET**.
2. Select [**DATE/TIME**] and press **CAL**.
3. Adjust the date (respectively year, month and day) by pressing **CAL** each time.
4. Adjust the time (respectively hour, minutes and seconds) by pressing **CAL** each time.

Recorder

1. Press **SET**.
2. Select [**4-20 mA**] and press **CAL**.
3. Select the desired range and press **CAL**.
4. Select the desired minimum level corresponding to a 4 mA recorder output and press **CAL**.
5. Select the desired maximum level corresponding to a 20 mA recorder output and press **CAL**.

Adjustment:

1. Select [**ADJUST**] and press **CAL**.
2. Connect a mA-meter to the 4-20 mA outputs.
3. Select **OUT1 / 4 mA**, press **CAL** and read the mA-meter.
4. Adjust to the correct current and press **CAL**.
5. Repeat point 3-4 for all other currents.

Password

1. Press **SET**.
2. Select [**PASSWORD**] and press **CAL**.
3. A private code can be programmed to avoid undesired access to the instrument. Enter your secret sequence of 5 keys and press **CAL**.

1. Press **SET**.
2. Select [CHANNELS] and press **CAL**.
3. Select [INPUT] and press **CAL**.
4. Select the desired range and press **CAL**.
5. Select [ATC] and press **CAL**.
6. Select the desired temperature compensation input and press **CAL**.
7. Select [CONTROL] and press **CAL**.
8. Select the control type (on/off or proportional).
9. Pre-set the desired low level (Lo).
10. Select the corresponding relay (REL).
11. Pre-set the desired high level (Hi).
12. Select the corresponding relay (REL).
13. Pre-set the desired hysteresis (Hy) between the on and off state of the relays.
14. Pre-set the desired delay (dt) after which the relays should be activated, when one of the pre-set levels is exceeded.

Proportional control only:

15. Pre-set the desired the proportional control area (dP), in which proportional regulation should start.
16. Pre-set the desired pulse-time (Pt) of the relays. The relays will pulse at a rate proportional to the regulation difference.

Optional alarm procedure:

17. Select [ALARM AFTER] and press **CAL**.
18. The display shows the alarm time after which the process control should be interrupted in case of a failure. Pre-set the desired alarm-time and press **CAL**.
19. Select [REL] and press **CAL**.
20. Select the desired relay to be closed in case of an alarm and press **CAL**.

Data-logging

1. Press **SET**.
2. Select [DATA LOGGER] and press **CAL**.
3. Select [INTERVAL] and press **CAL**.
4. Select the desired time interval between the data-logging and press **CAL**.

View the stored values:

1. Select [TABLE], [GRAPH] or [SEND] and press **CAL**.
2. Select the desired form to display the stored data or send them to a computer and press **CAL** to continue. The scales of a graph can be modified by pressing **CAL**. Follow the instructions on the screen.

Erase the stored values:

1. Select [ERASE] and press **CAL**. Follow the instructions on the screen.

RS485

1. Press **SET**.
2. Select [RS485] and press **CAL**.
3. Select [BAUDRATE] and press **CAL**.
4. Select the desired baudrate and press **CAL**.
5. Select [IDENTIF.NUMBER] and press **CAL**.
6. Select the desired identification number (0...999) and press **CAL**. The instrument can be identified by e.g. a computer when a specific number is allocated to it.

- A blinking decimal point warns you for unstable measurements. Wait to read the display!
- Stirring the solution during the measurements promotes the homogeneity and is obligatory!
- The instrument will refuse automatic calibration when the electrode is unstable. Insufficient stirring or a worn electrode may be the cause. See the *Technical Assistance Forms* at the end of this manual.

(P822)

Measurements

1. Should you want to recalibrate, select the desired measurement and press **CAL**. Follow the instructions on the screen. The displayed possibilities depend on the channel set-up and selected measurement.

pH calibration:

2. Select the proper buffer from one of the values in memory or enter manually a special value and press **CAL**. The unused buffers should be switched off.
3. Remove the electrodes from their fittings, rinse them with distilled water and immerse them in the first buffer solution. Select [CALIBRATE], press **CAL** and follow the instructions on the screen until the calibration is finished.
4. After rinsing the electrodes with distilled water, insert them back in their fittings and start the control.

mV calibration:

2. Remove the electrodes from their fittings, rinse them with distilled water, immerse them in a standard solution of known potential, adjust to the proper value and press **CAL**.
3. After rinsing the electrodes with distilled water, insert them back in their fittings and start the control.

Oxygen calibration:

2. The display shows the salinity correction. Select the proper values and press **CAL**. Leave salinity correction to zero unless you are going to measure in heavily salted solutions such as e.g. sea-water (35 g/l). Select [CALIBRATE], press **CAL** and follow the instructions on the screen until the calibration is finished.
3. Remove the electrode from its fitting and rinse it with distilled water. The electrode exposed to the air reaches an equilibrium corresponding to the partial pressure of oxygen and thus to saturation in water at the given temperature. The instrument shows the measured saturation, current, temperature, and will calibrate automatically when readings are stable.
4. After rinsing the electrode with distilled water, insert it back in its fitting and start the control. Stirring the solution during the measurements promotes the homogeneity and is **obligatory!** The advection rate must be at least 10 cm/s.

Temperature calibration:

1. Without Pt1000, adjust the manual temperature compensation and proceed by pressing **SET**. Should you want to recalibrate, press **CAL**.
2. Immerse the Pt1000 in a standard solution of known temperature. Calibrate to the proper value and press **CAL**.

Adjustment:

2. Measure the real value of the processed liquid with an other electrode and meter. Adjust to the real value and press **CAL**. Press **STOP** here if you wish to return to the original readings.
- Press **↑** to maximum measurement.
 - Press **↓** to minimum measurement.
 - Press **STOP** to stop the control without interrupting the measurements.

- A blinking decimal point warns you for unstable measurements. Wait to read the display!
- Stirring the solution during the measurements promotes the homogeneity and is obligatory!
- The instrument will refuse automatic calibration when the electrode is unstable. Insufficient stirring or a worn electrode may be the cause. See the *Technical Assistance Forms* at the end of this manual.

(P862)

Measurements

1. Should you want to recalibrate, select the desired measurement and press **CAL**. Follow the instructions on the screen. The displayed possibilities depend on the channel set-up and selected measurement.

Capacitive compensation:

2. Capacitive compensation increases accuracy in the very low conductivity ranges (<10 $\mu\text{S}/\text{cm}$). Select whether this compensation should be applied or not. Verify if the attached electrode is completely dry. Select [COMP.ON] or [COMP.OFF] and press **CAL**. Follow the instructions on the screen.

Conductivity calibration:

2. Select the proper standard from one of the values in memory or enter manually a special value and press **CAL**. Also select the temperature to which all future measurements will be referred to.
3. Remove the electrode from its fitting.
4. After rinsing the electrode several times with the standard solution, immerse it in that standard. The solution temperature is not so critical but should lie between 0°C and 30°C. Select [CALIBRATE], press **CAL** and follow the instructions on the screen until the calibration is finished.
5. After rinsing the electrode with distilled water, insert it back in its fitting and start the control.

Temperature calibration:

2. Without Pt1000, adjust the manual temperature compensation and proceed by pressing **SET**. Should you want to recalibrate, press **CAL**.
3. Immerse the Pt1000 in a standard solution of known temperature. Calibrate to the proper value and press **CAL**.

Adjustment:

2. Measure the real value of the processed liquid with an other electrode and meter. Adjust to the real value and press **CAL**. Press **STOP** here if you wish to return to the original readings.
- Press **↑** to show the maximum measurement.
 - Press **↓** to show the minimum measurement.
 - Press **STOP** to stop the control without interrupting the measurements.

A pH electrode is active and stable only after wetting! For this purpose it must be immersed for **at least ten hours** in a 3...4 M KCl solution. During short interruptions (e.g. storage) the electrode should be immersed in a 3...4 M KCl solution. In doing this it is always kept ready for use. When the interruption is longer than a month, refill the closing cap with 3...4 M KCl and plug it on the electrode tip in order to protect the glass bulb. Before use, ensure that the reference part of the electrode is topped up with a 3...4 M KCl solution.

Avoid a low pressure inside the electrode! Therefore always remove the closure from the refilling aperture during the measurements as well as during the calibration. This allows the saltbridge solution to flow through the ceramic liquid junction and prevents contamination of the electrolyte. For the same reason, the inside level should always be higher than the outside level of the measuring solution. Close the refilling aperture again when storing the electrode.

A polluted electrode may be cleaned with a soft detergent or 0.1 M HCl. Greasy substances may be removed with acetone or alcohol (**never do this with plastic electrodes!**).

If the electrode is polluted by proteinaceous materials (such as blood), it should stand in a cleaning solution overnight and then be cleaned with distilled water before use. The pH electrode wears away by being used. If the electrode tends to respond slower and calibration becomes difficult, even after cleaning, it should be replaced by a new one.

Maintenance of metal electrodes

Metal electrodes (Pt, Ag, Au): Metal electrodes are always ready for use. During short interruptions they are immersed in distilled water. **They should be cleaned regularly:**

- Silver electrodes are immersed in a concentrated ammonia solution during one hour.
- Platinum or gold electrodes are immersed in concentrated nitric acid during one hour.

Maintenance of conductivity electrodes

A conductivity cell is active and stable only after wetting! For this purpose it must be immersed for **at least one hour** in distilled water. Rinse the cell always after use and store it in distilled water (add some detergent to keep the spongy platinum surface in perfect condition).

A polluted cell may be cleaned with a soft detergent or diluted nitric acid. Greasy substances may be removed with acetone or alcohol (**never do this with plastic electrodes!**).

Maintenance of oxygen electrodes

A dissolved oxygen electrode is active and stable only when polarised! The electrode body has almost unlimited life and requires no maintenance. After prolonged use of the electrode, it may become deactivated. An indication is that the electrode no longer responds correctly to calibration. In this case:

1. Carefully screw off the measuring head.
 2. Cautiously remove the precipitated deposits with the aid of filter paper. **Do not use grinding paper or a glass fibre!**
 3. Rinse several times with distilled water and shake carefully off the water drops.
 4. Take the measuring head, fill it with new electrolyte solution and very slowly screw it onto the electrode while holding it vertically.
 5. **Let the electrode rest for a few hours!** The electrode is now ready for use again.
- *Replace the membrane only when damaged, not when calibration is no longer possible!*

1. Press **SET**.
2. Select [GLP] and press **CAL**.
3. Select [SHOW REPORT] and press **CAL**. Browse with **←→↑↓** to show a complete calibration report.
4. Select [SEND] and press **CAL** to send the report to a computer.

pH/mV/O₂/°C-controller P822

```
-----
Date           : 15/09/1998
Time          : 11:32:04
Version       : 1.0
```

SETTINGS

```
-----
Identification No : 003
Password         : OFF
GLP reminder     (h) : ---

Channel1        : pH
ATC             : °C1
Manual temp.    (°C) : 23.4

Channel2        : mV
```

pH1

CALIBRATION

```
-----
Date           : 15/09/1998
Time          : 11:20:06
```

Buffer 4.01/6.87

```
Slope          (%) : 98.4
Zero point     (pH) : 6.87
```

AVERAGE VALUES

```
-----
Response time  (s) : 11
Temperature    (°C) : 22.6
```

STATISTICS

```
-----
Slope          (%) : +0.4
Zero point     (pH) : +0.02
Response time  (s) : -1
```

mV2

```
Date           : 15/09/1998
Time          : 11:20:06
Calibr. at    (mV) : 235
```

Calibration reminder

1. Select [INTERVAL] and press **CAL**.
2. Select the desired time interval between each automatic warning for a new calibration of the electrodes and press **CAL**.

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| Specifications | | P822 |
|--------------------------|----------------------------|---|
| MODES | <i>pH</i> | 0...14 pH |
| | <i>mV</i> | ±2000 mV |
| | <i>Oxygen</i> | 0...60 mg/l (0...600%) |
| | <i>Temperature</i> | -30...+130 °C |
| RESOLUTION | <i>pH/mV</i> | 0.01 pH/ 1 mV |
| | <i>Oxygen</i> | 0.05 mg/l (0.5%) |
| | <i>Temperature</i> | 0.1 °C (1 °C: >99.9 °C) |
| CHANNELS | <i>pH/mV/O₂</i> | 2 |
| | <i>Temperature</i> | 2 for Pt1000 (or Pt100) |
| INPUTS | | screw terminals |
| TEMPERATURE COMPENSATION | | -30...+130 °C (O ₂ : 0...50 °C) |
| SALINITY COMPENSATION | <i>Oxygen</i> | 0...40 ppt |
| ACCURACY | <i>pH/mV</i> | 0.1% ± 1 digit |
| | <i>Oxygen</i> | 1% ± 1 digit |
| CALIBRATION | <i>Temperature</i> | 0.3 °C |
| | <i>Reminder</i> | 0...999 h |
| | <i>pH</i> | 1...2 points |
| | <i>Buffers</i> | 9 pre-programmed 2 user specified |
| | <i>ISO-pH</i> | 6...8 pH |
| | <i>Slope</i> | 80...120% |
| | <i>mV/Oxygen</i> | 1 point |
| | <i>Temperature</i> | ±10 °C |
| | <i>On/Off</i> | ✓ |
| CONTROL | <i>Proportional</i> | ✓ |
| | <i>LCD</i> | 128x64 pixels |
| DISPLAY | | 128x64 pixels |
| GLP PROTOCOL | | ✓ |
| REAL TIME CLOCK | | ✓ |
| RECORDER | <i>Built-in</i> | 1200 points |
| | <i>Interval</i> | 1...9999 s |
| | <i>Two Outputs</i> | 4...20 mA, programmable max. 300 Ω load |
| RS485 | <i>Baud Rate</i> | 300...9600 b/s |
| | <i>Computer</i> | up to 31 units via RS485/RS232 interface |
| DATA-LOGGING | <i>Values</i> | 1200 + °C/date/time |
| | <i>Interval</i> | 1...9999 s |
| RELAY | | 4 solid state |
| OUTPUTS | | 12...250 VAC/ 0.5 A |
| HELP FUNCTION | <i>Languages</i> | EN, NL, FR, DE |
| IDENTIFICATION NUMBER | | ✓ |
| PASSWORD PROTECTION | | ✓ |
| POWER SUPPLY | <i>Mains</i> | 210...250 VAC |
| DIMENSIONS | <i>WxDxH</i> | 20x9x12 cm |
| WEIGHT | | 1.2 kg |

| Specifications | | P862 |
|-----------------------|-----------------------|---|
| MODES | <i>Conductivity</i> | 0...2000 mS/cm |
| | <i>Resistivity</i> | 0...20 MΩ.cm |
| | <i>TDS</i> | 0...100 g/l |
| | <i>Temperature</i> | -30...+130 °C |
| RESOLUTION | <i>Conductivity</i> | 0.01 µS/cm |
| | <i>Resistivity</i> | 1 Ω.cm |
| | <i>TDS</i> | 0.01 mg/l |
| | <i>Temperature</i> | 0.1 °C |
| CHANNELS | <i>Conductivity</i> | 2 |
| | <i>Temperature</i> | 2 for Pt1000 (or Pt100) |
| INPUTS | | screw terminals |
| TEMPERATURE COMP. | | -30...+130 °C, aut. or man. |
| ACCURACY | <i>Conductivity</i> | 0.5% f.s. of range |
| | <i>Temperature</i> | 0.3 °C |
| CALIBRATION | <i>Reminder</i> | 0...999 h |
| | <i>Conductivity</i> | 1 point |
| | <i>Cell Constant</i> | 0.1/1/10 cm ⁻¹ ± 50% |
| | <i>Temp. Coef.</i> | natural waters (EN27888) |
| | <i>Ref. Temp.</i> | 20° or 25 °C |
| | <i>Cal. Solutions</i> | 0.01/0.1/1 M KCl |
| | <i>Cap. Comp.</i> | ✓ |
| | <i>Temperature</i> | ±10 °C |
| CONTROL | <i>On/Off</i> | ✓ |
| | <i>Proportional</i> | ✓ |
| DISPLAY | <i>LCD</i> | 128x64 pixels |
| GLP PROTOCOL | | ✓ |
| REAL TIME CLOCK | | ✓ |
| RECORDER | <i>Built-in</i> | 1200 points |
| | <i>Interval</i> | 1...9999 s |
| | <i>Two Outputs</i> | 4...20 mA, programmable max. 300 Ω load |
| RS485 | <i>Baud Rate</i> | 300...9600 b/s |
| | <i>Computer</i> | up to 31 units via RS485/RS232 interface |
| DATA-LOGGING | <i>Values</i> | 1200 + °C/date/time |
| | <i>Interval</i> | 1...9999 s |
| RELAY | | 3 solid state |
| OUTPUTS | | 12...250 VAC/ 0.5 A |
| HELP FUNCTION | <i>Languages</i> | EN, NL, FR, DE |
| IDENTIFICATION NUMBER | | ✓ |
| PASSWORD PROTECTION | | ✓ |
| POWER SUPPLY | <i>Mains</i> | 210...250 VAC |
| DIMENSIONS | <i>WxDxH</i> | 20x9x12 cm |
| WEIGHT | | 1.2 kg |

| Temp. | 1.68 | 4.01 | 6.87 | 9.18 | 12.45 | <i>NIST buffers</i> |
|--------------|-----------------|-----------------|-----------------|-----------------|------------------|---------------------|
| 0 °C | 1.666 pH | 4.003 pH | 6.984 pH | 9.464 pH | 13.423 pH | |
| 5 °C | 1.668 pH | 3.999 pH | 6.951 pH | 9.395 pH | 13.207 pH | |
| 10 °C | 1.670 pH | 3.998 pH | 6.923 pH | 9.332 pH | 13.003 pH | |
| 15 °C | 1.672 pH | 3.999 pH | 6.900 pH | 9.276 pH | 12.810 pH | |
| 20 °C | 1.675 pH | 4.002 pH | 6.881 pH | 9.225 pH | 12.627 pH | |
| 25 °C | 1.679 pH | 4.008 pH | 6.865 pH | 9.180 pH | 12.454 pH | |
| 30 °C | 1.683 pH | 4.015 pH | 6.853 pH | 9.139 pH | 12.289 pH | |
| 40 °C | 1.694 pH | 4.035 pH | 6.838 pH | 9.068 pH | 11.984 pH | |
| 50 °C | 1.707 pH | 4.060 pH | 6.833 pH | 9.011 pH | 11.705 pH | |
| 60 °C | 1.723 pH | 4.091 pH | 6.836 pH | 8.962 pH | 11.449 pH | |
| 70 °C | 1.743 pH | 4.126 pH | 6.845 pH | 8.921 pH | | |
| 80 °C | 1.766 pH | 4.164 pH | 6.859 pH | 8.885 pH | | |

| Temp. | 4.00 | 7.00 | 9.21 | 10.01 | <i>Standard buffers</i> |
|--------------|----------------|----------------|----------------|-----------------|-------------------------|
| 5 °C | 3.99 pH | 7.08 pH | 9.45 pH | 10.24 pH | |
| 10 °C | 3.99 pH | 7.06 pH | 9.38 pH | 10.18 pH | |
| 20 °C | 3.99 pH | 7.02 pH | 9.26 pH | 10.06 pH | |
| 25 °C | 4.00 pH | 7.00 pH | 9.21 pH | 10.01 pH | |
| 30 °C | 4.01 pH | 6.99 pH | 9.16 pH | 9.87 pH | |
| 40 °C | 4.03 pH | 6.98 pH | 9.06 pH | 9.89 pH | |
| 50 °C | 4.06 pH | 6.97 pH | 8.99 pH | 9.83 pH | |
| 60 °C | 4.09 pH | 6.98 pH | 8.93 pH | 9.79 pH | |

| Temp. | 0.01 M | 0.1 M | 1 M | <i>Conductivity standards</i> |
|--------------|----------------|-----------------|-----------------|-------------------------------|
| 0 °C | 776 µS | 7.15 mS | 65.4 mS | |
| 5 °C | 896 µS | 8.22 mS | 74.1 mS | |
| 10 °C | 1020 µS | 9.33 mS | 83.2 mS | |
| 15 °C | 1147 µS | 10.48 mS | 92.5 mS | |
| 20 °C | 1278 µS | 11.67 mS | 102.1 mS | |
| 25 °C | 1413 µS | 12.88 mS | 111.8 mS | |
| 30 °C | 1552 µS | 14.12 mS | 121.7 mS | |

| Temp. | 3 M KCl | <i>Redox standard</i> |
|-------|---------|-----------------------|
| 0 °C | 224 mV | |
| 5 °C | 219 mV | |
| 10 °C | 214 mV | |
| 15 °C | 212 mV | |
| 20 °C | 204 mV | |
| 25 °C | 199 mV | |
| 30 °C | 194 mV | |
| 35 °C | 189 mV | |
| 40 °C | 184 mV | |
| 45 °C | 179 mV | |
| 50 °C | 174 mV | |