

# DISSOLUTION ACCESSORIES

## *Dissolution Media Concentrates*

***Dissolution Accessories offers a complete line of Dissolution Media Concentrates.***

***Diluting the concentrate by simply adding purified water will create “ready to use” dissolution media within a few minutes.***



**DILUTE DISSOLUTION MEDIA CONCENTRATE – MIX – RUN DISSOLUTION TEST**

Our pre-mixed concentrates significantly reduce the amount of time and cost to prepare dissolution media. The preparation of dissolution media is a time consuming procedure on a daily basis; by using the dissolution media concentrates you will reduce the preparation time by up to 75%, so laboratory analysts/chemists can focus on performing more valued added activities.

All dissolution media concentrates are available in convenient sizes, which can be used to produce 6L, 10L or 25L dissolution media. This makes it the ideal solution for pharmaceutical research, formulation and quality control laboratories.

The small bottle sizes will reduce storage space requirements and lowers freight charges.

DA dissolution media concentrates are available in different bottle sizes (250ml, 500ml and 1L) of supplying the right solution for individual dissolution testing needs.

**The following DA dissolution media concentrates are available:**

Product	250 ml Bottle Size (230.8 ml Concentrate Volume) Dilutes to 6L of media	250 ml Bottle Size (250,0 ml Concentrate Volume) Dilutes to 10L of media	1 L Bottle Size (961.5 ml Concentrate Volume) Dilutes to 25L of media
0.01 N HCl	DMC020-06	DMC020-10	DMC020-25
0.1 N HCl	DMC021-06	DMC021-10	DMC021-25
Acetate Buffer, pH 4.5	DMC022-06	DMC022-10	DMC022-25
Potassium Phosphate, pH 5.8	DMC023-06	DMC023-10	DMC023-25
Potassium Phosphate, pH 6.0	DMC024-06	DMC024-10	DMC024-25
Potassium Phosphate, pH 6.8	DMC025-06	DMC025-10	DMC025-25
Potassium Phosphate, pH 7.2	DMC026-06	DMC026-10	DMC026-25
Potassium Phosphate, pH 7.4	DMC027-06	DMC027-10	DMC027-25
Potassium Phosphate, pH 7.5	DMC028-06	DMC028-10	DMC028-25
SLS 0.50% *	DMC029-06	DMC029-10	
Simulated Gastric Fluid (without enzyme)	DMC030-06	DMC030-10	DMC030-25
Simulated Intestinal Fluid (without enzyme)	DMC031-06	DMC031-10	DMC031-25

\* SLS 0.50% 500ml bottle size  
(400ml of concentrate volume produces 6L of diluted media and  
500ml of concentrate volume produces 10L of diluted media)



Please visit our webshop for the complete product range of Dissolution Accessories.

Order online at: [www.dissolutionaccessories.com](http://www.dissolutionaccessories.com)

Dissolution Accessories • PO Box 173 • 4900 AD Oosterhout • The Netherlands

Tel: +31(0)162-471485 • Fax: +31(0)162-471486

[info@dissolutionaccessories.com](mailto:info@dissolutionaccessories.com)