

## SENSELINE

WORKING IN THE FIELD BECOMES A PLEASURE

## MODEL SENSELINE F410

## pH/mV/°C Meter



- Clear display of measured values and settings
- Watertight and dustproof to IP67 specifications
- Continue indication of electrode and battery status
- Automatic endpoint and buffer recognition
- Automatic Temperature Compensation
- Data memory for 30 GLP data points
- 3-point calibration with predefined or user-defined buffers

**TOTALLY WATERPROOF  
BY DESIGN**

## SPECIFICATIONS

F410	Measuring range	Resolution	Accuracy
pH	0.00 ~ 14.00	0.01	± 0.01
mV	-1999 ~ 1999	1	± 1
Temperature	-5.0 ~ 105°C	0.1°C	± 0.5°C
Sensor inputs	BNC (>10 <sup>12</sup> Ω); NTC 30 KΩ (both IP67)		
Power requirements	4AA batteries 1.5V or NiMH accumulators 1.3V		
Operation environment	0 ~ 40°C, 5 ~ 85% rel. humidity (non-cond.)		
Size / Weight	220 x 90 x 45mm / 325g (without batteries)		

## ORDERING INFORMATION

<b>F410</b>	Meter only option includes F410 pH portable meter, operation manual and four AA batteries
<b>F410T</b>	Meter kit option includes F410 + 2 pH buffers + KCl solution + epoxy pH electrode with built-in ATC (IP67) and carrying case